




Arthrex Angel™ System

## Bone Marrow Concentration

Indication-specific PRP and PRF gel preparations from bone marrow aspirate



ARTHREX  
**angel**  
SYSTEM


**Arthrex**®

# Bone Marrow Concentration

Bone marrow contains nucleated and progenitor cells such as hematopoietic and mesenchymal stem cells. The Arthrex® Angel BMC Kit (ABS-10072) for Arthrex Angel™ System is a convenient and rapid means of concentrating the cellular contents and growth factors in bone marrow aspirate: \*

	Platelet Concentration (K/ $\mu$ l)	Nucleated Cell Concentration (K/ $\mu$ l)	Hematopoietic Cell Concentration (K/ $\mu$ l)	Total Neutrophil ( $\times 10^6$ )
<b>BMA</b>	87.7 $\pm$ 6.4	24.5 $\pm$ 15.6	0.002 $\pm$ 0.001	612.1
<b>BMC</b>	787.0 $\pm$ 317.6	240.5 $\pm$ 186.6	0.081 $\pm$ 0.056	132.9
<b>Increase Above Baseline</b>	~9x	~10x	~33x	↓80%

\* These values are averages, which may differ in some cases. Starting volume was always 60 ml anticoagulated BMA from iliac crest, HCT 7%. The average volume was 1.8 ml of BMC.




CFUs after 96h of incubation time

## Harvesting



### Calcaneus Harvest Technique

Make a small incision 1 cm anterior and 1 cm plantar to the insertion of the Achilles tendon over the lateral portion of the calcaneus, taking care to avoid the sural nerve. When inserting the needle do not exceed a depth of 3 cm. Aspirate a small volume of bone marrow redirecting as necessary until the desired volume of BMA is obtained.



### Posterior Iliac Crest Harvest Technique

Make a small incision at the desired location over the PSIS of the iliac crest. Use the needle tip to locate the center of the iliac crest. Insert the needle and advance 3 cm. Aspirate the bone marrow slowly, redirecting as necessary. Repeat until the desired volume is obtained.



### Arthroscopic Distal Femur Harvest Technique

Bone marrow aspiration should occur before drilling tunnels. Arthroscopically insert the needle in the apex of the femoral notch to a depth of 3 cm. Turn off the arthroscopic fluid before removing the trocar and attaching the syringe. Slowly aspirate the bone marrow. In order to obtain the desired volume it may be necessary to rotate 90° or withdraw the needle 0.5 cm when aspirating; prevent withdrawing the needle past the 2 cm mark.





### Arthroscopic Proximal Humerus Harvest Technique

Bone marrow aspiration should occur before any fixation implants are inserted. Arthroscopically insert the needle in the location where the first anchor would be placed; do not exceed a depth of 3 cm. Turn off the arthroscopic fluid before removing the trocar and attaching the syringe. Aspirate the bone marrow slowly. In order to obtain the desired volume it may be necessary to rotate 90° or withdraw the needle 0.5 cm when aspirating; prevent withdrawing the needle past the 2 cm mark.

## Applications

Concentrated bone marrow aspirate is a bioactive fluid that supports the ingrowth and resorption of a variety of biocompatible and resorbable scaffolds.<sup>1,2,3,4</sup> Autograft and allograft but also synthetic material can be mixed with concentrated bone marrow aspirate prior to application to an orthopaedic site:





Arthrex® Angel BMC-Kit

- 1 Platelet Cuvette
- 2 Rotating Valve
- 3 Variable Volume Separation Chamber
- 4 Three-Compartment Reservoir Bag
- 5 60ml Wrapped Specimen Cup
- 6 Male-Female Luer Lock Plugs
- 7 Whole Blood Spike Adapter
- 8 20ml Luer Lock Syringe
- 9 Syringe-Activated PRP Valve
- 10 BMA Clot Filter
- 11 Fenestrated Bone Marrow Needle
- 12 Vac Lok Syringe, 30ml
- 13 Hypodermic Needle, 20G

## Ordering Information

Arthrex Angel™-System  
 Arthrex® Angel BMC Kit with Bone Marrow Needle  
 ACD-A (20ml)

ABS-10066  
 ABS-10072  
 SAAV222.G00

### References:

1. Orii H, Sotome S, Chen J, et al. Beta-tricalcium phosphate (beta-TCP) graft combined with bone marrow stromal cells (MSCs) for posterolateral spine fusion. *J Med Dent Sci* 2005;52:51-7.
2. Kai T, Shao-qing G, Geng-ting D. In vivo evaluation of bone marrow stromal-derived osteoblasts-porous calcium phosphate ceramic composites as bone graft substitute for lumbar intervertebral spinal fusion. *Spine* 2003;28:1653-8.
3. Betsch M, Schneppendahl J, Thuns S, et al. Bone Marrow Aspiration Concentrate and Platelet Rich Plasma for Osteochondral Repair in a Porcine Osteochondral Defect Model. *PLoS ONE* 2013; 8(8): e71602.
4. Shih H, Shih L, Sung T, Chang Y. Restoration of bone defect and enhancement of bone ingrowth using partially demineralized bone matrix and marrow stromal cells. *Journal of Orthopaedic Research* 2005; 23(6): 1293-9.



www.arthrex.com

***This brochure is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. The contents of this brochure do not constitute a recommendation for specific treatment. As part of this professional usage, the medical professionals must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professionals should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's Directions For Use. Within the framework of medical freedom and based on the risk-benefit analysis for the patient in question, the treating physician bears sole responsibility for the appropriateness of treating the patient with the Arthrex Angel™ System.***